

Preparation of stool samples

Swiss TPH

Swiss Tropical and Public Health Institute Schweizerisches Tropen- und Public Health-Institut Institut Tropical et de Santé Publique Suisse






3/4



Pack and close in the plastic cover provided.

- Fill the empty container ¾ full with stool (approx. the size of an apricot)
- Seal tightly
- · Label with the date of collection










Send it in the prestamped, preaddressed bags together with the completed order form.

Stool samples in SAF fixative solution (Container with transparent liquid)







Pack and close in the plastic sleeve provided.

- A CAUTION: Do not touch or drink the liquid! It is hazardous to your health!
- Add a hazelnut-sized amount of stool to the liquid using the wooden spatula and mix well
- Seal tightly
- Label with the date of collection

If multiple tests are required, please collect the samples on different days.

The samples can be sent separately or together.

Information about the tests: www.swisstph.ch