

Swiss Tropical and Public Health Institute Schweizerisches Tropen- und Public Health-Institut Institut Tropical et de Santé Publique Suisse

Associated Institute of the University of Basel


Ingrid Felger
Molecular Diagnostics Unit
Dept. of Medical Parasitology & Infection Biology

Diagnostics for malaria surveillance

- Parasite detection
- Quantification
- Genotyping



Programmatic decisions rely on surveillance data



Standardized tools & metrics

understand infectious reservoir for malariameasure changes in malaria transmission



Malaria infections

Symptomatic



requires diagnostic sensitivity of 200 parasites/µL blood



Microscopy and some RDTs


good LM & good RDT:

≈ 100 parasites/µL blood

asymptomatic → submicroscopic


Relevant for transmission?

All to be found?
All to be treated?
Proportional contribution to transmission?
Changes with levels of endemicity?



Malaria infections

Symptomatic → asymptomatic → submicroscopic



on average 50%

70-80% in low transmission20% in high transmission

P. vivax problematic!

5-8 x lower mean densities by qPCR & LM compared to Pf in asymptomatic individuals

Bousema et al. Nature Rev. Microbiol 2014

The asymptomatic Reservoir of Infection



Parasite prevalence

- simple field work
- cross-sectional surveys
- Correlates well with exposure /_{mol}FOI



optimize test sensitivity to explore submicroscopic reservoir




Development of <u>ultra-sensitive DNA-based</u> qPCR

- > TARE-2: telomere-associated repetitive element 2
- ➤ 1.6 kb long blocks of 10-12 **135-bp repeat units** with slightly degenerate sequences
- > approx. 250–280 copies / genome
- var-ATS: acidic terminal segment (semi-conserved)
- > 59 *var* genes in 3D7 genome


Good correlation between TARE2, ATS and standard 18S rRNA qPCR

→ useful for quantification of low densities





Implication for prevalence rate: plus 16%



Major advantage of multicopy PCR target:

Pooling of 10 blood samples without loss of sensitivity

P. falciparum prevalencein 498 individuals from Tanzania

Hofmann et al. PLOS Med 2015

The asymptomatic Reservoir of Infection




Parasite prevalence


Gametocyte prevalence

- ➤ Detected on transcript level by qRT-PCR or NASBA (Pett et al. 2016 Mal J)
- Most sensitive molec. marker: pfs25 & pvs25
- ➤ RNA based → difficult sampling




Implication for prevalence rate: plus 16%






P. falciparum prevalence in 498 individuals from Tanzania

Proportion of gametocyte carriers by pfs25 qRT-PCR





Parasite prevalence

Gametocyte density / prevalence

Infectivity to mosquitoes

Membrane feeding assays

More data needed

- Pf and Pv
- Different endemic sites

Landmark paper

Churcher et al. 2013 Elife



Do we need molecular assays for gametocytes?

Detection?

No

- ☐ All infections produce gametocytes at a steady rate (Eksi et al.2012 PLOS Pathogen)
- □ Detection is a matter of blood volume

Quantification?


Yes - for certain purposes

- □ Research: Describing the asymptomatic infective reservoir
- Monitoring of drug effects in clinical trials in comparator groups
- Vaccine effects on gametocytes in comparator groups



Current developments: increase sensitivity of assays

- ultra-sensitive qPCR (PNG, Brazil, Thailand)
 - high copy number targets for Pf and Pv
 - □ qRT-PCR targeting 18S ribosomal RNA



- increased blood volume 3 mL venous bleed (in PNG; 4 malaria species)
 - □ white blood cell depletion → concentration of PCR template
 - □ magnetic enrichment of gametocytes
 - □ male / female gametocyte quantification
- evaluation of iso-thermal methods / suitable for point-of-care tools
 - □ LAMP (Loop-mediated isothermal **AMP**lification)
 - □ NASBA (Nucleic Acid Sequence Based Amplification) RNA-based

"Paper Machine" for POC Molecular Diagnostics



Disposable device (5x10cm) intended for point-of-care use in resource-limited environments

Prototype tested 5 cells *E.coli* in 50 µL human serum

Device that integrates sample preparation and **loop-mediated isothermal amplification (LAMP; 1hr @ 65°C)** with end point detection using a hand-held UV source and camera phone



Elimination phase: Detection of malaria infections

Surveillance – response strategy

1 - 3 - 7

Day 1:

Diagnosis
Index case
Treatment
Digital
Reporting

Day 3: Epidemiol. Evaluation at CDC Day 7: Response in foci Screen



Focal Screen and Treat or MSAT


Example China Elimination until 2020

Weizhong Yang 2015 ECTMIH Conference




Molec. Barcodes
MS haplotype
Deep Seq Data















SNP-based genotyping methods

Barcoding


genome-wide single-nucleotide-polymorphisms are amplified by PCR

- ~20 SNPs, each detected by qPCR
- Barcodes for regional or global clone discrimination



Chromosome

1 *P. falciparum* genome = **14** chromosomes




Alternative SNP-based genotyping methods

Barcoding


genome-wide single-nucleotide-polymorphisms are amplified by PCR

- ~20 SNPs each detected by qPCR
- Work load!
- Barcodes for regional or global clone discrimination



Problem to reconstruct haplotype in case of multi-clone infections


► used primarily for **single clone** infections/low endemicity only



Amplicon "deep sequencing"


A. Lerch unpublished

This blood sample 1 harbours 2 different parasite clones (= 2 genomes)




PCR reaction 1

Multiplex PCR with several sequence-tagged primer pairs



Amplicon "deep sequencing"



Example shows 2/384 blood samples with 5 markers/genome sequenced



Conclusion

human infectious reservoir







Infectiousness to mosquitoes

- ✓ Well described
- ✓ Low-density reservoir

Not enough knowledge on contribution of low-density infections

Needed: More hard work → MFAs

Needed

- Diagnostics Development for surveillance tools / POC
- Spatial and temporal heterogeneity of the infectious reservoir



Nov. 2016 Madang, PNG

Lincoln Timinao informing school children about Membrane Feeding Assays



Information of school teachers at the IMR Entomology laboratory, Madang, PNG



IMR field team and Natalie Hofmann screening PNG schoolchildren November 2016

Thanks





Natalie Hofmann
Rahel Wampfler
Felistas Mwingira
Sarah Javati
Maria Gruenberg
Anita Lerch
Clara Moniz
Tom Smith



Ivo Mueller Cristian Koepfli



Leanne Robinson Lincoln Timinao Stephan Karl



Seif Shekalaghe

Financial support

Swiss National Science Foundation ICEMR (Internatl. Centers of Excellence in Malaria Research) Bill and Melinda Gates Foundation

