

Mobile Phone and Handheld Microscopy for Helminth Diagnoses

Swiss TPH Winter Symposium – Dec. 7, 2017

Isaac I. Bogoch, MD, MS, FRCPC, DTM&H
Division of Infectious Diseases, Toronto General Hospital
Assistant Professor, University of Toronto, Canada

Outline

1. Why use **mobile phone** and handheld microscopes?
2. What **devices** have been used to date?
3. What are some of the **limitations**?
4. What are the **future** directions to overcome limitations?



Thank you

- **Dr. Jean Coulibaly** – Université Félix Houphouët-Boigny, Ivory Coast
- **Mr. Shaali Ame** – Public Health Laboratory, Pemba, Tanzania
- **Prof. Juerg Utzinger** – Swiss Tropical and Public Health Institute
- **Prof. Jennifer Keizer** – Swiss Tropical and Public Health Institute
- **Prof. Aydogan Ozcan** – University of California, Los Angeles
- **Prof. Dan Fletcher** – University of California, Berkeley
- **Dr. Richard Kobina Dadze** – University of Cape Coast, Ghana
- **Dr. Jason Andrews** – Stanford University

- **Funding:** Grand Challenges Canada



Princeton, ~2003

Laos, 2014



***“Praziquantel
Cheers”***

I. Why use mobile phone and handheld microscopes?



The Issues

- Lack of **diagnostic support** at clinical and public health levels in resource-constrained settings
 - Especially in rural areas
- Some facilities can check **basic** blood tests
- Smaller places and more rural clinics have **few or no** diagnostic tests
 - Send out specimens vs. transfer people to larger clinics



So what happens?

- Transferring people and/or specimens:
 - Time
 - Money
 - Sick people
 - No time - too sick, or have to be home to help support family
 - No money











- We attempt to deliver **equitable healthcare** and **public health** by bringing **quality laboratory diagnostics** to low-resource settings
- ***Bring the lab to the people rather than moving patients or specimens***
- Goal: develop, validate, and implement **portable, easy-to-use** microscopes to areas most in need
- Converting **mobile phones** into microscopes and using other **handheld microscope devices**

Why do you want a mobile phone?

- Portable, **easy** to use, **battery** powered
- Digitize an image
 - Send to to others via **text** or **email**
 - Automated diagnoses with **computer vision** and **machine learning**
- Record other **health information**/data with image
- Mark the **GPS** coordinates


Such technology could be used for rapid, portable, point-of-contact diagnoses for clinical or public health applications in resource-constrained settings




2. What devices have been used to date?



A long time ago in a galaxy far,
far away....







T. trichiura



A. lumbricoides








Terrible!


Anton van Leeuwenhoek, 1600s



Foldscope



Reversed Lens Microscope



Device	Organism Detected	Sample Size	Sn/Sp (%)	NPV/PPV (%)
Ball lens mounted to mobile phone	<i>A. lumbricoides</i> , <i>T. trichiura</i> , hookworm	199	69.4/61.5	92.3/23.2
Foldscope	<i>S. haematobium</i>	49	55.9/93.3	95.0/48.3
Reversed-lens Cellscope	<i>S. haematobium</i>	49	67.6/100.0	100.0/57.7



Handheld and Mobile Phone Microscopes


- Handheld microscope with mobile phone attachment
- Weighs 0.5 kg
- Good success field studies on STH and schistosomiasis



Newton Nm¹










Device	Organism Detected	Sample Size	Sn/Sp (%)	NPV/PPV (%)
Ball lens mounted to mobile phone	<i>A. lumbricoides</i> , <i>T. trichiura</i> , hookworm	199	69.4/61.5	92.3/23.2
Foldscope	<i>S. haematobium</i>	49	55.9/93.3	95.0/48.3
Reversed-lens Cellscope	<i>S. haematobium</i>	49	67.6/100.0	100.0/57.7
Newton Nm I with mobile phone attached	<i>S. mansoni</i>	226	91.7/99.5	91.7/99.5
Newton Nm I with mobile phone attached	<i>S. haematobium</i>	226	81.1/97.1	94.8/88.6
3D printed mobile phone microscope	<i>S. haematobium</i>	59	72.1/100.0	100.0/57.1

S. haematobium: sensitivity by egg count



S. haematobium
<50 low intensity
≥50 high intensity



- ***How do these devices work in real-world setting by people who actually use them in their day-to-day clinical or public health work?***
- We have evaluated the performance of devices, used by local microscopists, integrated in community-based screening and treatment programs for malaria, schistosomiasis and soil-transmitted helminth infection















The unintended “selfie”



WHEN MICROSCOPE SLIDES



**DESCRIBE YOUR
RELATIONSHIP STATUS**



3. What are some of the limitations?



What are some of the limitations?

- Design of devices:
 - **Optics** are good in most devices
 - Better design devices for **high-throughput** of samples
 - may need to re-think sample preparation
- Implementation of devices:
 - Need validation in **real-world** settings
 - Issues with scale-up when devices are “*good enough*”

4. What are the future directions to overcome limitations?



Capabilities prototype devices?

- Quality imaging at the sub-micron level
- Computer vision and machine learning
 - Automated identification and quantification of organisms
- Tracking GPS coordinates to map diagnoses
- Ability to input other meaningful clinical data
- Endless possibility for other applications



*A word on cost-effectiveness

The cost of the device is trivial as long as it is robust, has sufficient optics, and is designed for efficient throughput




Summary

- Mobile Phone microscopy has incredible **potential** to improve the **quality of care** in *clinical* and *public health* settings
- Validate new devices in **real-world** settings
- **Scale up** training and usage in appropriate settings

References

- Bogoch II, Andrews JR, Speich B, Utzinger J, Ame SM, Ali SM, Keiser J. Mobile phone microscopy for the diagnosis of soil-transmitted helminth infections: a proof-of-concept study. *Am J Trop Med Hyg*. 2013 Apr;88(4):626-9.
- Linder E, Grote A, Varjo S, Linder N, Lebbad M, Lundin M, Diwan V, Hannuksela J, Lundin J. On-chip imaging of *Schistosoma haematobium* eggs in urine for diagnosis by computer vision. *PLoS Negl Trop Dis*. 2013 Dec 5;7(12):e2547.
- Koydemir HC, Gorocs Z, Tseng D, Cortazar B, Feng S, Chan RY, Burbano J, McLeod E, Ozcan A. Rapid imaging, detection and quantification of *Giardia lamblia* cysts using mobile-phone based fluorescent microscopy and machine learning. *Lab Chip*. 2015 Feb 16;15(5):1284-93.
- Wei Q, Luo W, Chiang S, Kappel T, Mejia C, Tseng D, Chan RY, Yan E, Qi H, Shabbir F, Ozkan H, Feng S, Ozcan A. Imaging and sizing of single DNA molecules on a mobile phone. *ACS Nano*. 2014 Dec 23;8(12):12725-33.
- Bogoch II, Coulibaly JT, Andrews JR, Speich B, Keiser J, Stothard JR, N'goran EK, Utzinger J. Evaluation of portable microscopic devices for the diagnosis of *Schistosoma* and soil-transmitted helminth infection. *Parasitology*. 2014 Dec;141(14):1811-8.
- Cybulski JS, Clements J, Prakash M. Foldscope: origami-based paper microscope. *PLoS One*. 2014 Jun 18;9(6):e98781.
- Switz NA, D'Ambrosio MV, Fletcher DA. Low-cost mobile phone microscopy with a reversed mobile phone camera lens. *PLoS One*. 2014 May 22;9(5):e95330.
- WHO (2002). Prevention and control of schistosomiasis and soil-transmitted helminthiasis: a report of a WHO expert committee. WHO Technical Report Series 912, 1–57.





Thank you

isaac.bogoch@uhn.ca



Costs

- Build robust microscopes that **function** well
- **Resist** the urge to build the **cheapest** microscopes that sacrifice image quality
- We have seen costs for typical components drop and many of these items are available in developing countries
 - Smart phones
 - Webcams



Costs – an exercise

- Lifespan of microscope = 10 years
- Cost of laboratory technician
 - \$4000/year
 - x10 years
 - = **\$40,000**
- Works
 - 48 weeks per year
 - 5 days per week



Costs – an exercise

▪ \$1 device

- Processes **30** samples per day
- 7,200 samples per year
- 72,000 over 10 years
- Total cost over 10 years
 $\$40,000 + \$1 = \$40,001$

▪ **Cost per sample: \$0.56**

▪ \$500 device

- Process **40** samples per day
- 9,600 samples per year
- 96,000 samples over 10 years
- Total cost over 10 years
 $\$40,000 + \$500 = \$40,500$

▪ **Total cost per sample: \$0.42**



Costs – an exercise

- Image quality with \$500 device likely better than \$1 device
- Added 'soft' costs with cheaper device - lowered sensitivity and specificity
 - Missed diagnoses, greater healthcare expenses
 - Time off of work, morbidity, etc.
- ***The cost of the device is trivial, throughput and proper diagnosis is everything***